



# Wireless Veterinary Ultrasound

C10EQ

Rectal probe for  
reproductive diagnostic  
in cows and mares

**KONTED**

[www.konted.cn](http://www.konted.cn)






## Characteristics:

- Wifi connection between probe and your phone  
No cable, easier to carry and use.
- Whole probe is waterproof, convenient for cleaning
- The probe is combined with a auxiliary rod.

## Specifications:

- Monitor: Smart phone/tablet PC Series
- Working System: Android/IOS
- Frequency: 6.5-8 MHz
- Scan mode: Linear
- Elements: 128
- Grey levels: 256 levels
- Scan depth: 40/60/80/100 mm
- Display mode: B mode
- Pseudo-color: 8 mode  
(Gray coding/Red saturation/Green saturation/Blue saturation/Yellow saturation/Purple saturation/Rainbow coding/ Metal coding)
- Software: Gain adjustment, Depth, 8 TCG, Enhancement, Dynamic Range, Cine-Loop, Grid, save images and video
- Wi-Fi : 802.11n/2.4G/5G
- Battery: 6800mAh
- Working time: 6 hours (14 hours in stand by)
- Main unit size: 780x30x40 mm
- Weight: 690 g



128 Elements Rectal probe

Length:78cm

Built-in Wi-Fi and Lithium battery